

[www.thundertilt.com](http://www.thundertilt.com)

**Platinum Power Standing Wheelchair**  
**Fits Most Sunrise, Invacare and Pride Mid Wheel Drive Bases**  
**Also Custom Fit to Most Centre Wheel Drive and Rear Wheel Drive Bases**



**Seating Drive Position  
Shown on Sunrise QM 710**



**Standing Drive Position  
Shown on Sunrise QM 710**



**Power Tilt Option  
Shown on Sunrise QM 710**



**Power Recline with Shear Reduction  
Shown on Invacare TDX SP**

## Options shown on Quantum 6000 Base



**Standing Drive position**

### Features Include on All Standing Systems:

- Power Chair Can Be Driven In the Standing Position if the Chair has the Required Electronics
- Adjustable Shear Reduction for Travel of Back and Footplate
- Adjustable Depth Knee Block
- Right and Left Knee Block Independently Adjustable Angle
- Adjustable Chest Belt
- Swing Back Adjustable Height Armrests Desk or Full Length
- Adjustable Speed for Power Stand, Power Tilt and Power Recline if the Chair has the Required Electronics

### Options & Accessories

- Platinum Power Tilt
- Platinum Power Recline with Shear Reduction
- Platinum Power Leg rest
- Swing Away Joystick Mount
- Lock on Joystick Knob
- R-Net Electronic Available for all Power Bases
- Throw Joystick Operation
- Attachment Points

**Platinum Power Standing Wheelchair**  
**Custom Fit to Most Rear Wheel Drive Bases**

**Options Shown on a Storm Base**



**Additional Function & Operation of the Rear Wheel Standing System**

- Front Castor are Locked in Forward Position for Stability
- Centre Mount Leg Rest Powers Down to Allow Standing
- Centre Mount Leg Rest Powers up to Clear Castor for Driving
- Power Leg Rest will also Elevate to 60 Degree Power Elevating Leg Rest Function
- To Operate Leg and or Stand Functions, Front Castor must be in Lock Position